




## ecoSHIELD\* Natural Nitrile PF 250



- Category III PPE glove (PPE Directive 89/686/EEC)
- Complex Design - For mortal and irreversible risks
- Powder-free translucent nitrile glove
- Ambidextrous
- 250 mm / 0.08 mm (EN 420:2003 + A1:2009)
- AQL 1.5 (EN 374-2:2014 Level 2)
- Viral penetration test (ISO 16604:2004 Procedure B & ASTM F1671-97b)
- Waterproof and for low chemical protection
- Tested for chemical permeation (EN 16523-1:2015 supersedes EN 374-3:2003)
- Type I hypersensitivity eliminated - Type IV hypersensitivity reduced

## ecoSHIELD\* Eco Nitrile PF 250



- Category III PPE glove (PPE Directive 89/686/EEC)
- Complex Design - For mortal and irreversible risks
- Powder-free green/white nitrile glove
- twinSHIELD\* double-walled protection
- Ambidextrous
- 250 mm / 0.10 mm (XL/10: 260 mm) (EN 420:2003 + A1:2009)
- AQL 0.65 (EN 374-2:2014 Level 3)
- Viral penetration test (ISO 16604:2004 Procedure B & ASTM F1671-97b)
- Waterproof and for low chemical protection
- Tested for chemical permeation (EN 16523-1:2015 supersedes EN 374-3:2003)
- Type I hypersensitivity eliminated - Type IV hypersensitivity reduced

### ecoSHIELD\* Natural Nitrile PF 250

5329-14-6 Sulfamic Acid, 15%	Not Tested
64742-49-0 Petroleum benzene 80-100°C	Not Tested
66-81-9 Cycloheximide	Not Tested

### ecoSHIELD\* Eco Nitrile PF 250

5329-14-6 Sulfamic Acid, 15%	<b>LEVEL 6</b> 480 min
64742-49-0 Petroleum benzene 80-100°C	<b>LEVEL 6</b> 480 min
66-81-9 Cycloheximide	<b>LEVEL 6</b> 480 min

## ecoSHIELD\* Natural Nitrile PF 250

Mixed Solution Hydranal® -Composite 2	Not Tested
Mixed Solution Mucocit®-T branded mixture	Not Tested
10127-02-3 Acridine orange	Not Tested
77-86-1 Tris (hydroxymethyl) aminomethane Sat. solution	Not Tested
127-09-3 Sodium acetate Sat. solution	Not Tested
7664-93-9 Sulphuric Acid 50%	LEVEL 6 480 min
108-95-2 Phenol 0.1% solution	Not Tested
75-59-2 Tetramethylammonium hydroxide 2.5%	Not Tested
7664-93-9 Sulphuric Acid 10%	Not Tested
7681-52-9 Sodium Hypochlorite 13%	Not Tested
1310-58-3 Potassium Hydroxide 40%	Not Tested
1310-73-2 Sodium Hydroxide, 50%	Not Tested
1310-73-2 Sodium Hydroxide 40%	LEVEL 6 480 min
7664-38-2 Phosphoric Acid, 30%	Not Tested
7722-84-1 Hydrogen Peroxide 30%	Not Tested

## ecoSHIELD\* Eco Nitrile PF 250

Mixed Solution Hydranal® -Composite 2	LEVEL 6 480 min
Mixed Solution Mucocit®-T branded mixture	LEVEL 6 480 min
10127-02-3 Acridine orange	LEVEL 6 480 min
77-86-1 Tris (hydroxymethyl) aminomethane Sat. solution	LEVEL 6 480 min
127-09-3 Sodium acetate Sat. solution	LEVEL 6 480 min
7664-93-9 Sulphuric Acid 50%	LEVEL 6 480 min
108-95-2 Phenol 0.1% solution	LEVEL 6 480 min
75-59-2 Tetramethylammonium hydroxide 2.5%	LEVEL 6 480 min
7664-93-9 Sulphuric Acid 10%	LEVEL 6 480 min
7681-52-9 Sodium Hypochlorite 13%	LEVEL 6 480 min
1310-58-3 Potassium Hydroxide 40%	LEVEL 6 480 min
1310-73-2 Sodium Hydroxide, 50%	LEVEL 6 480 min
1310-73-2 Sodium Hydroxide 40%	LEVEL 6 480 min DR 10%
7664-38-2 Phosphoric Acid, 30%	LEVEL 6 480 min
7722-84-1 Hydrogen Peroxide 30%	LEVEL 6 480 min

## ecoSHIELD\* Natural Nitrile PF 250

111-30-8 Glutaraldehyde 25%	LEVEL 6 480 min
50-00-0 Formaldehyde 37%	Not Tested
1239-45-8 Ethidium Bromide 5%	LEVEL 6 480 min
7722-64-7 Potassium permanganate 6.5%	Not Tested
7803-57-8 Hydrazine monohydrate 98%	Not Tested
7647-01-0 Hydrochloric Acid 37%	LEVEL 1 22 min
67-56-1 Methanol, 20%	Not Tested
124-18-5 n-Decane	Not Tested
67-63-0 Isopropanol 70%	LEVEL 2 32 min
108-87-2 Methylcyclohexane 99,9%	Not Tested
Mixed Solution Bacillol 30 Foam	LEVEL 0 0 min
110-82-7 Cyclohexane	Not Tested
79-06-1 Acrylamide 40%	LEVEL 6 480 min
142-82-5 n-Heptane 99%	Not Tested
78-83-1 Isobutanol 99%	Not Tested

## ecoSHIELD\* Eco Nitrile PF 250

111-30-8 Glutaraldehyde 25%	Not Tested
50-00-0 Formaldehyde 37%	LEVEL 6 480 min DR 25%
1239-45-8 Ethidium Bromide 5%	LEVEL 6 480 min
7722-64-7 Potassium permanganate 6.5%	LEVEL 4 153 min
7803-57-8 Hydrazine monohydrate 98%	LEVEL 4 150 min
7647-01-0 Hydrochloric Acid 37%	LEVEL 3 98 min
67-56-1 Methanol, 20%	LEVEL 3 65 min
124-18-5 n-Decane	LEVEL 2 58 min
67-63-0 Isopropanol 70%	LEVEL 2 58 min
108-87-2 Methylcyclohexane 99,9%	LEVEL 2 55 min
Mixed Solution Bacillol 30 Foam	LEVEL 2 48 min
110-82-7 Cyclohexane	LEVEL 6 480 min
79-06-1 Acrylamide 40%	LEVEL 6 480 min
142-82-5 n-Heptane 99%	LEVEL 2 41 min DR 48%
78-83-1 Isobutanol 99%	LEVEL 2 39 min

**ecoSHIELD\* Natural Nitrile PF 250**

67-63-0 Isopropanol 100%	<b>LEVEL 1</b> 21 min
64-17-5 Ethanol 70%	<b>LEVEL 1</b> 11 min
598-75-4 Secondary isoamyl alcohol 98%	Not Tested
110-54-3 n-Hexane 95%	Not Tested
71-36-3 Butanol 99,4%	Not Tested
Mixed Solution Bacillol AF	Not Tested
108-95-2 Phenol 50%	Not Tested
Mixed Solution Terralin protect	Not Tested
75-12-7 Formamide 99%	Not Tested
62-53-3 Aniline 99,9%	Not Tested
97-88-1 Butyl methacrylate 99,9%	Not Tested
67-68-5 Dimethyl Sulfoxide 99% (DMSO)	Not Tested
64-17-5 Ethanol 99.8%	<b>LEVEL 0</b> 7 min
7664-39-3 Hydrofluoric Acid 40%	Not Tested
64-19-7 Acetic Acid 100%	<b>LEVEL 0</b> 3 min

**ecoSHIELD\* Eco Nitrile PF 250**

67-63-0 Isopropanol 100%	<b>LEVEL 2</b> 33 min
64-17-5 Ethanol 70%	<b>LEVEL 2</b> 31 min
598-75-4 Secondary isoamyl alcohol 98%	<b>LEVEL 2</b> 30 min
110-54-3 n-Hexane 95%	<b>LEVEL 1</b> 28 min
71-36-3 Butanol 99,4%	<b>LEVEL 1</b> 25 min
Mixed Solution Bacillol AF	<b>LEVEL 1</b> 22 min
108-95-2 Phenol 50%	<b>LEVEL 1</b> 15 min
Mixed Solution Terralin protect	<b>LEVEL 1</b> 13 min
75-12-7 Formamide 99%	<b>LEVEL 1</b> 11 min
62-53-3 Aniline 99,9%	<b>LEVEL 1</b> 10 min
97-88-1 Butyl methacrylate 99,9%	<b>LEVEL 1</b> 10 min
67-68-5 Dimethyl Sulfoxide 99% (DMSO)	<b>LEVEL 1</b> 10 min
64-17-5 Ethanol 99.8%	<b>LEVEL 0</b> 9 min
7664-39-3 Hydrofluoric Acid 40%	<b>LEVEL 0</b> 8 min
64-19-7 Acetic Acid 100%	<b>LEVEL 0</b> 4 min

## ecoSHIELD\* Natural Nitrile PF 250

7697-37-2 Nitric Acid 70%	LEVEL 0 2 min
64-18-6 Formic acid 97.8%	Not Tested
1330-20-7 Xylene 98,5%	LEVEL 0 0 min
1336-21-6 Ammonium Hydroxide 25%	Not Tested
60-24-2 2-Mercaptoethanol 99%	Not Tested
68-12-2 Dimethyl Formamide 99%	Not Tested
584-84-9 Toluene diisocyanate 95%	Not Tested
75-56-9 Propylene oxide 99%	Not Tested
7664-93-9 Sulphuric Acid 95%-98%	Not Tested
60-29-7 Diethyl ether 99%	Not Tested
75-15-0 Carbon Disulfide 99,9%	Not Tested
76-05-1 Trifluoroacetic acid 99%	Not Tested
67-64-1 Acetone 99,8%	LEVEL 0 0 min
109-99-9 Tetrahydrofuran 99.9%	Not Tested
100-42-5 Styrene 99.9%	Not Tested

## ecoSHIELD\* Eco Nitrile PF 250

7697-37-2 Nitric Acid 70%	LEVEL 0 4 min
64-18-6 Formic acid 97.8%	LEVEL 0 3 min
1330-20-7 Xylene 98,5%	LEVEL 0 2 min
1336-21-6 Ammonium Hydroxide 25%	LEVEL 1 16 min
60-24-2 2-Mercaptoethanol 99%	LEVEL 0 1 min
68-12-2 Dimethyl Formamide 99%	LEVEL 0 1 min
584-84-9 Toluene diisocyanate 95%	LEVEL 0 0 min
75-56-9 Propylene oxide 99%	LEVEL 0 0 min
7664-93-9 Sulphuric Acid 95%-98%	LEVEL 0 0 min
60-29-7 Diethyl ether 99%	LEVEL 0 0 min
75-15-0 Carbon Disulfide 99,9%	LEVEL 0 0 min
76-05-1 Trifluoroacetic acid 99%	LEVEL 0 0 min
67-64-1 Acetone 99,8%	LEVEL 0 0 min
109-99-9 Tetrahydrofuran 99.9%	LEVEL 0 0 min
100-42-5 Styrene 99.9%	LEVEL 0 0 min

## ecoSHIELD\* Natural Nitrile PF 250

110-86-1 Pyridine	Not Tested
75-09-2 Dichloromethane 99%	LEVEL 0 0 min
67-56-1 Methanol 99,9%	LEVEL 0 0 min
7664-39-3 Hydrofluoric Acid 48%	Not Tested
64-18-6 Formic acid 98,5%	Not Tested
141-78-6 Ethyl Acetate 99.8%	LEVEL 0 0 min
109-89-7 Diethylamine 99,5%	LEVEL 0 0 min
67-66-3 Chloroform 99,8%	LEVEL 0 0 min
107-13-1 Acrylonitrile 99%	Not Tested
75-05-8 Acetonitrile 99,9%	LEVEL 0 0 min

## ecoSHIELD\* Eco Nitrile PF 250

110-86-1 Pyridine	LEVEL 0 0 min
75-09-2 Dichloromethane 99%	LEVEL 0 0 min
67-56-1 Methanol 99,9%	LEVEL 0 0 min
7664-39-3 Hydrofluoric Acid 48%	LEVEL 0 0 min
64-18-6 Formic acid 98,5%	LEVEL 0 0 min
141-78-6 Ethyl Acetate 99.8%	LEVEL 0 0 min
109-89-7 Diethylamine 99,5%	LEVEL 0 0 min
67-66-3 Chloroform 99,8%	LEVEL 0 0 min
107-13-1 Acrylonitrile 99%	LEVEL 0 0 min
75-05-8 Acetonitrile 99,9%	LEVEL 0 0 min